





# Zoll R Series Testing and Readiness

The 4P communication tool	
<p><b>Purpose:</b> <i>Why are we doing this?</i></p>	<p>Team members must know that the Zoll R Series defibrillator is ready for use after successful completion of the most recent code readiness test (manual or automatic) as indicated by a green check in the code readiness display window.</p>
<p><b>Picture:</b> <i>What will it look like when we get there?</i></p>	<p>The code readiness display window located on the front panel of the Zoll R Series should have a green check mark indicating the defibrillator is ready for patient use. If there is a red X, this indicates something triggered a failure during code readiness testing, and the process <u>must</u> be repeated. The red X indicates the defibrillator is <u>not ready</u> for patient use and may have malfunctioned. If the red X persists after further testing (see below), the defibrillator must be removed from use and a ticket should be submitted to Clinical Engineering for service.</p> <div style="display: flex; justify-content: space-around;">   </div>
<p><b>Plan:</b> <i>What is our plan to change?</i></p>	<p>All team members will be educated regarding the meaning of the <span style="color: green;">✓</span> and the <span style="color: red;">✗</span>.</p>
<p><b>Part:</b> <i>What is my role going to be in the plan?</i></p>	<p>Designated team members will continue to manually test the defibrillator (as specified by each unit) and ensure it is operating properly as indicated by the green check.</p> <p>All team members are responsible for ensuring a defibrillator is ready for use, prior to patient use, by visualizing the green check in the code readiness display window.</p> <p>If there is a red X, switch the defibrillator to Monitor mode&gt; select Report Data&gt; select Test Log&gt; view and correct the cause of the failure as indicated (the failure cause will be highlighted). If the red X persists, follow the process above and engage Clinical Engineering. <u>Note:</u> R Series Plus users may need to put the monitor in manual mode to troubleshoot.</p> <p>If an error message appears, follow the directions provided to correct and re-test. If unable to resolve, remove from use and engage Clinical Engineering.</p> <p>Contact your department leader with questions.</p>