







DISINFECTION & FOOD SERVICE WIPES

	Product	Kill/Wet/ Contact Time *	Recommended Surfaces	Category of Disinfectant	Dilution expiration	Lawson #	Notes
	Sani-Cloth AF3 (Gray top)	3 minutes	Hard surfaces - IV poles, BP cuffs, wheel chairs, commodes. Must be used to clean Abbott glucometers.	Quaternary Ammonium. Product is alcohol free	See manufacturer's expiration on container (set at 2 years)	# 99403	For general disinfection; not for use with isolation precautions. Wear gloves.
	Sani-Cloth Bleach (Orange top)	4 minutes	Compatible with hard surfaces ** Specific equipment like CareFusion Alaris pumps.	Sodium hypochlorite (Bleach)	See manufacturer's expiration on container (set at 11 month expiration)	# 99404	For all isolations, including <i>C.diff</i> , <i>Norovirus</i> . Wear gloves.
	Sani-wipe (Red top)	1 minute	Food service items -dining tables, coffee pots, refrigerators, microwaves, office desks/mouse/keyboard, toys.	Quaternary Ammonium. Product is alcohol free	See manufacturer's expiration on container	# 592232	Gloves not needed.
	PeriSept	2 minutes (Hepatitis B/C in 5 minutes)	Not recommended for copper, brass, granite, marble, or zinc. Nor used on unsealed/uncoated marble or terrazzo floors.	Peroxyacetic acid / hydrogen peroxide	24 hours	Used only by EVS	General use disinfectant for EVS.
	NorOxyCdiff	2 minutes (Hepatitis B/C in 5 minutes)	Not recommended for copper, brass, granite, marble, or zinc. Nor used on unsealed/uncoated marble or terrazzo floors.	Peroxyacetic acid / hydrogen peroxide	Test before use	Used only by EVS	Used by EVS for 'misting'

*Surface must stay wet for the specified time. Do not wipe dry.

**To potentially decrease risk of cracking vinyl, wipe surfaces after kill time with water and a clean cloth. Wear gloves.

DISINFECTION & FOOD SERVICE WIPES

	Product	Kill/Wet/ Contact Time *	Recommended Surfaces	Category of Disinfectant	Dilution expiration	Lawson #	Notes
	Virex II 256	10 minutes	Compatible with most medical device material and hard surfaces (Check with manufacturer of equipment)	Quaternary ammonium	24 hrs	Used by EVS	Used by EVS for disinfecting surfaces in patients rooms

*Surface must stay wet for the specified time. Do not wipe dry.

**To potentially decrease risk of cracking vinyl, wipe surfaces after kill time with water and a clean cloth. Wear gloves.